


 GeneSafe™

The first non-invasive prenatal test that
screens for single-gene disorders

GeneSafe™ the evolution of NIPT

A non-invasive prenatal test that screens multiple genes for mutations causing severe genetic disorders in the fetus



GeneSafe™


analyses circulating cell-free fetal DNA (cfDNA) from a maternal blood sample.

The test is performed after **10 weeks** of pregnancy.

GeneSafe™ works as a **complementary screen** to traditional and genome-wide NIPT **PrenatalSafe® KARYO**. It screens for several life-altering genetic disorders that are not screened with current NIPT technology, allowing a **complete picture** of the risk of a pregnancy being affected by a genetic disorder.



facilitates early diagnosis of single-gene disorders.



It involves 3 different levels of screening:



This test screens for **5 common inherited recessive genetic disorders**, such as **Cystic Fibrosis, Beta-Thalassemia, Sickle cell anaemia, Deafness autosomal recessive type 1A, Deafness autosomal recessive type 1B.**

Genes screened: CFTR, CX26 (GJB2), CX30 (GJB6), HBB



This test screens for **44 severe genetic disorders** due to **de novo** mutations (a gene mutation that is not inherited) in **25 genes**

Genes screened: ASXL1, BRAF, CBL, CHD7, COL1A1, COL1A2, COL2A1, FGFR2, FGFR3, HDAC8, JAG1, KRAS, MAP2K1, MAP2K2, MECP2, NIPBL, NRAS, NSD1, PTPN11, RAF1, RIT1, SETBP1, SHOC2, SIX3, SOS1



This test screens for both **inherited** and **de novo** single-gene disorders and represents a combination of the tests  and  providing a complete picture of the pregnancy risk.



allows detection of common inherited genetic disorders in the fetus

GENE	GENETIC DISORDER
CFTR	Cystic Fibrosis
CX26 (GJB2)	Deafness autosomal recessive type 1A
CX30 (GJB6)	Deafness autosomal recessive type 1B
HBB	Beta-Thalassemia
HBB	Sickle cell anemia

The inherited recessive disorders screened by  are the most common in the European population



identifies fetal conditions that could be missed by traditional prenatal screening.

GENE	SYNDROMIC DISORDERS
JAG1	Alagille syndrome
CHD7	CHARGE syndrome
HDAC8	Cornelia de Lange syndrome 5
NIPBL	Cornelia de Lange syndrome 1
MECP2	Rett syndrome
NSD1	Sotos syndrome 1
ASXL1	Bohring-Opitz syndrome
SETBP1	Schinzel-Giedion syndrome
SIX3	Holoprosencephaly
NOONAN SYNDROMES	
BRAF	Cardiofaciocutaneous syndrome 1
CBL	Noonan syndrome-like disorder with or without juvenile myelomonocytic leukemia (NSLL)
KRAS	Noonan syndrome/cancers
MAP2K1	Cardiofaciocutaneous syndrome 3
MAP2K2	Cardiofaciocutaneous syndrome 4
NRAS	Noonan syndrome 6/cancers
PTPN11	Noonan syndrome 1/ LEOPARD syndrome/cancers
PTPN11	Juvenile myelomonocytic leukemia (JMML)
RAF1	Noonan syndrome 5/LEOPARD syndrome 2
RIT1	Noonan syndrome 8
SHOC2	Noonan syndrome-like disorder with loose anagen hair
SOS1	Noonan syndrome 4


GENE	SKELETAL DISORDERS	
COL2A1	Achondrogenesis, type II or hypochondrogenesis	
	Achondroplasia	
	CATSHL syndrome	
FGFR3	Crouzon syndrome with acanthosis nigricans	
	Hypochondroplasia	
	Muenke syndrome	
	Thanatophoric dysplasia, type I	
COL1A1	Thanatophoric dysplasia, type II	
	Ehlers-Danlos syndrome, classic	
	Ehlers-Danlos syndrome, type VIIA	
	Osteogenesis imperfecta, type I	
COL1A2	Osteogenesis imperfecta, type II	
	Osteogenesis imperfecta, type III	
	Osteogenesis imperfecta, type IV	
	Ehlers-Danlos syndrome, cardiac valvular form	
COL1A2	Ehlers-Danlos syndrome, type VIIB	
	Osteogenesis imperfecta, type II	
	Osteogenesis imperfecta, type III	
COL1A2	Osteogenesis imperfecta, type IV	
	CRANIOSYNOSTOSIS SYNDROMES	
	FGFR2	Antley-Bixler syndrome without genital anomalies or disordered steroidogenesis
Apert syndrome		
Crouzon syndrome		
Jackson-Weiss syndrome		
Pfeiffer syndrome type 1		
Pfeiffer syndrome type 2		
Pfeiffer syndrome type 3		

GeneSafe™ detects de novo mutations in 25 genes causing 44 different genetic disorders. The genetic conditions screened by this innovative test are often present in the fetus in the absence of a family history of the condition. This is a paradigm shift in prenatal screening. GeneSafe™ screens for de novo mutations that cannot be detected by standard carrier screening, as these mutations are not present in parents' somatic cells. The genetic disorders screened by GeneSafe™ can cause **skeletal dysplasias, cardiac defects¹⁻³, multiple congenital anomalies⁴⁻⁵, autism⁶, epilepsy⁷, and/or intellectual disability⁸⁻⁹.**

1. Homsy J, et al. *Science*. 2015;350:1262–6. 2. Zaidi S, et al. *Nature*. 2013;498:220–3. 3. Sifrim A, et al. *Nat Genet*. 2016;48:1060–5. 4. Ng SB, et al. *Nat Genet*. 2010;42:790–3. 5. Hoischen A, et al. *Nat Genet*. 2011;43:729–31. 6. O’Roak BJ, et al. *Nat Genet*. 2011;43:585–9. 7. Allen AS, et al. *Nature*. 2013;501:217–21. 8. de Ligt J, et al. *N Engl J Med*. 2012;367:1921–9. 9. Rauch A, et al. *Lancet*. 2012;380:1674–82.



screens for conditions common across all maternal ages



All pregnant women – **regardless of age** – are at equal risk of the genetic conditions screened by GeneSafe™. Although the occurrence of each disorder is relatively rare, the cumulative rate of occurrence of these conditions (**~1 in 600**) is similar to that of Down Syndrome, in younger women¹⁰.

10. McRae J, et al. Prevalence and architecture of de novo mutations in developmental disorders. *Nature* 2017, 542, 433–438



can identify conditions that may have otherwise gone undetected until after birth or into childhood



Many disorders screened with  are not typically associated with abnormal prenatal ultrasound findings (especially in the first trimester), or may not be evident until late second/ third trimester, when confirmatory invasive testing can pose a risk of preterm birth, or baby's health after delivery.




screens for genetic disorders associated with advanced paternal age



While traditional NIPT screens for conditions typically associated with advanced maternal age (e.g. Down Syndrome), GeneSafe™ DE NOVO screens for genetic disorders (e.g. **Achondroplasia, Pfeiffer syndrome, Crozen syndrome, Apert syndrome, Thanatophoric dysplasia, Osteogenesis Imperfecta**, etc.), which are associated with **advanced paternal age** (men who are **>40 years old**)¹¹ ensuring a comprehensive screen for couple of advanced age. These disorders typically are caused by **mutations** arising during **spermatogenesis**. As a man ages, the chance for these errors to occur substantially increases.

¹¹. Kong A, et al.: Rate of de novo mutations and the importance of father's age to disease risk. *Nature* 2012, 488:471-475.

GeneSafe™ clear test results reporting



POSITIVE

Pathogenic / Likely Pathogenic mutation(s) detected

This result shows that the test detected one or more mutations in one or more genes.

A patient with a positive  GeneSafe™ test result should be referred for genetic counselling and should always be followed-up with an invasive diagnostic test for confirmation of test results, before any medical decisions are made.

Only **known pathogenic** and **likely pathogenic** mutations are reported.



NEGATIVE


NO Pathogenic / Likely Pathogenic mutation(s) detected

This result shows the test has not detected any disease causing mutation in the targeted genes screened.


Negative screening results mean that there is a **very low risk** that the fetus has one of the disorders screened with  GeneSafe™ although no guarantee may be given that the fetus is actually healthy.

 GeneSafe™ a groundbreaking technology allowing for a genetic analysis that is revolutionary

cfDNA isolation
from maternal plasma



Fetal Fraction
determination



Next Generation
Sequencing (NGS)

cfDNA mutation analysis

Advanced bioinformatic analysis to identify de novo mutations



 GeneSafe™ **Superior accuracy**



High Resolution: coverage >550X

HD

Exceptional performance:
Sensitivity and specificity >99%



Accurate measurement of Fetal Fraction (FF)



Low limit of detection: highly accurate at **low cfDNA quantity** (FF > 2%)



Low incidence of inconclusive results (<1%)



GeneSafe™ Test Characteristics

SIMPLE: a simple **blood sample (8-10 ml)** collected at **10+ weeks** of gestation is required



SAFE: it is a **non-invasive** test, **no risk** to the fetus and the mother



RELIABLE:
Sensitivity and specificity >99%



FAST:
Turnaround time of **10 days**



GeneSafe™ **Indication for testing**



Is intended for patients who meet any of the following criteria:

- *Advanced paternal age (men who are > 40 years old)*
- *Abnormal ultrasound finding(s) suggestive of monogenic disorder*
- *Patients wishing to avoid an invasive diagnostic procedure*
- *Patients at risk for genetic conditions screened*

The test is suitable for:

- *both single and twin pregnancies.*
- *patients whose pregnancies have been achieved by IVF techniques, including pregnancies with egg donation or surrogacy.*

GeneSafe™ 5 easy steps



Order the **test**



Fill in and Sign the **Test Requisition and Consent Form.**



Provide a **maternal blood sample (10 ml)**
at **10+ weeks** of pregnancy




Call **DHL** on **0844 248 0844** to arrange
collection of the sample, **that is it!**



Receive results in just **10 days**

GeneSafe™ Complimentary services



GeneSafe™ includes

Free **follow-up** of abnormal results

Free CVS or Amniocentesis
in collaboration with reference gynecologists

Reimbursement of the test fee
for cases with inconclusive test results

Genetic Counselling for Positive Results



advanced molecular diagnostics solutions using state-of-the art technologies



Test performed in Italy
(Rome or Milan)



Fast TAT: 10 days



20 years experience in prenatal
molecular diagnostics



Personalised genetic counselling
with genetic counselors experts
in discussing genetic test results
and familial risks.



Fully **Automated work-flow**, from
cfDNA extraction to data analysis



Test available **worldwide**



Over **200,000** genetic tests/year



Dedicated R&D team
Numerous peer-reviewed papers
published in renowned international
journals



Eurofins Biomnis UK Ltd., 17 Doman Road, Camberley, GU15 3DF
Laboratori e Studi Medici, Via Castel Giubileo, 11 - 00138 Roma
Laboratori e Studi Medici, Via Enrico Cialdini, 16 - 20161 Milano Tel.:
+39 06 8811270 | Fax: +39 0664492025

www.laboratoriogenoma.eu
www.prenatalsafe.co.uk
www.prenatalsafekaryo.co.uk
www.genesafe.co.uk

Next Step

Please call or email to find out more
UK toll free number: 0808 1691022
Email: NIPT@eurofins.co.uk



Genoma