

## Abstract/Introduction

Pharmaceutical manufacturing equipment must be properly cleaned to ensure the removal of product residue, cleaning chemical residue and microbes prior to manufacturing. Cleaning methods need to be developed and validated to prevent the risk of cross contaminated and adulteration of products. These methods need to be validated to confirm that the defined cleaning process sufficiently removes potential contaminants. Establishment of method limits and selection of the appropriate cleaning techniques and detection methods are critical to prove that the defined method conditions effectively clean the manufacturing equipment. Validation of the defined method conditions provides confidence in the defined method conditions. This poster will discuss the selection of the appropriate cleaning procedure, including selection of various sampling techniques, and the detection options available to monitor the amount of residual contaminant. The use of the appropriate tests during method validation are critical in proving that the method performs as intended, and these validation tests are discussed as well. In addition, revalidation of analytical techniques and the use of correction factors will be discussed.

## Establishment of Appropriate Limits

In the US FDA Guide to Inspections Validation of Cleaning Processes, it is stated, "The firm's rational for the residue limits established should be logical based on the manufacturer's knowledge of the materials involved and be practical, achievable, and verifiable." The FDA gives examples of analytical detection levels such as 10 ppm, biological activity levels such as 1/1000 (0.1%) of the normal therapeutic dose and organoleptic levels such as no visible residue. One can easily apply these examples to determine the amount of allowable carry-over of product residues.

- No more than 10 ppm of any product should appear in another product.
- No more than 0.1% of the normal therapeutic dose of any product will appear in the maximum daily dose of the following product.
- No amount of product residue should be visible on the surface of the equipment after the cleaning procedure has been performed.

## DISCUSSION

Visual Determination of Residue		Methods of Detection: There are multiple detection options for cleaning validation/verification samples producing quantitative information on residues.		Summary of Typical Validation Components*			
<p>The FDA states, "If the cleaning process is used only between batches of the same product (or different lots of the same intermediate in a bulk process), the firm need only meet a criteria of 'visibly clean' for the equipment."</p> <p>Enhance visual inspection by:</p> <ul style="list-style-type: none"><li>• Spiking coupons with different known amounts of residue</li><li>• Training personnel observe coupons to determine at which level coupon appears clean</li><li>• Acceptance level set at highest level of visual coupon that appears clean!</li></ul> <p><b>Disadvantages to Visual Inspection</b></p> <ul style="list-style-type: none"><li>• Too many variables that can influence results</li><li>• Coupons must be observed in same viewing conditions as equipment in the field</li><li>• Not all equipment can be viewed similar to coupon sitting on lab bench</li><li>• The lighting must be the same</li><li>• The angle of viewing must be the same</li><li>• The direction the viewer is from the surface must be the same</li><li>• Observer may not be able to see edge of stain instead of the body of the stain itself!</li><li>• Results are not quantitative</li></ul>		<p><b>Instrumentation</b></p> <p><b>IMS – Ion Mobility Spectrometry</b></p> <ul style="list-style-type: none"><li>• Characterizes chemical substances based on their gas phase ion mobilities</li><li>• Provides detection and quantitation of trace analytes</li><li>• Utilizes atmospheric pressure chemical ionization (APCI) – a soft ionization technique that produces molecular weight information</li></ul> <p><b>TOC – Total Organic Carbon</b></p> <ul style="list-style-type: none"><li>• Analysis is specific to organic compounds</li><li>• Theoretically measures all the covalently bonded carbon in water!</li></ul> <p><b>UV-Visible Spectrophotometry</b></p> <ul style="list-style-type: none"><li>• Commonly used for detection of small molecule active pharmaceutical ingredients or detergent residues for swab and rinse samples</li></ul> <p><b>HPLC - High Performance Liquid Chromatography</b></p> <ul style="list-style-type: none"><li>• Used for detection of small molecule active pharmaceutical ingredients or detergent residues for both swab and rinse samples allowing for separation of multiple components</li></ul> <p><b>GC and GC/MS – Gas Chromatography and Mass Spectrometry</b></p> <ul style="list-style-type: none"><li>• Used mainly for detection of detergent residue</li><li>• Specific to volatile and semi-volatile organic compounds</li></ul> <p><b>Swabbing</b></p> <ul style="list-style-type: none"><li>• Ideal for residues not easily removed with water rinsing</li><li>• Swabs can physically remove insoluble residues</li><li>• Chosen swabs must exhibit the following characteristics:<ul style="list-style-type: none"><li>1) Ability to recover desired residue from given surface</li><li>2) Ability to release the residue into an extraction solution for analysis</li><li>3) Must not contribute excessive interference or background against analysis</li></ul></li></ul>		<p><b>Benefits</b></p> <ul style="list-style-type: none"><li>• Ultra-fast quantitative analysis (~30 seconds per sample)</li><li>• Sub-nanogram sensitivity</li><li>• The ability to analyze a broad range of compounds with no chromatophore needed</li><li>• No mobile phases, columns or vacuum is required for operation</li><li>• Sample introduction via either thermal desorption or by high performance injection</li></ul> <p><b>Drawbacks</b></p> <ul style="list-style-type: none"><li>• Compounds must be water-soluble and ionizable for IMS detection to be used!</li><li>• Samples must be relatively clean</li><li>• Ultra-pure extraction solutions should be used</li><li>• Not suitable for multiple component matrices</li><li>• TOC analysis incorporates all organic molecules in solution</li><li>• Material may comprise carbon from various components and not just just compound of interest!</li><li>• Contaminating materials! need to be organic and contain carbon that can be oxidized under TOC test conditions!</li><li>• Samples must be water soluble (above MCL)</li><li>• Sensitive to interferences</li><li>• Lacks peak separation</li><li>• Chromophore required for specificity</li></ul> <p><b>Validation Component</b></p> <p><b>Accuracy</b></p> <p>Accuracy should be assessed at a minimum of 3 concentration levels, each prepared in triplicate. Typically performed with concentrations ranging from 80 to 120% of final theoretical concentration. Accuracy may include evaluation of recover from spiked swabs or from spiked surfaces.</p> <p><b>Precision</b></p> <p>Precision is frequently performed in conjunction with Accuracy. For precision the 100% spike level (whether swabs or surfaces) is prepared as 6 replicates. These 6 replicates are prepared by 2 separate analysts and the results are compared.</p> <p><b>Linearity</b></p> <p>Performed using a minimum of 5 concentration levels, each injected in duplicate. Linearity may range from the limit of quantitation up to 200% of the MCL.</p> <p><b>Specificity</b></p> <p>Typically both swabs and surfaces are evaluated to determine if interferences with the compound of interest are present.</p> <p><b>LOD</b></p> <p>Prepare standard solutions at the estimated LOD (3 preparations) and analyze.</p> <p><b>LOQ</b></p> <p>Prepare standard solutions at the estimated LOQ (3 preparations injected in duplicate) and analyze.</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>The ability of a standard or sample preparation solution to meet method specifications over time.</p>		<p><b>"Typical" Acceptance Criteria</b></p> <ul style="list-style-type: none"><li>• Swab recovery: mean recovery of 90% - 110% theoretical, %RSD <math>\leq</math> 10%</li><li>• Surface recovery: mean recovery of 85% - 115% theoretical, %RSD <math>\leq</math> 15%</li><li>• All system suitability meets method criteria.</li><li>• Intermediate precision data and precision data (combined) must have an RSD of <math>\leq</math> 15%.</li></ul> <p><b>Linearity</b></p> <p>Correlation coefficient (r) <math>\geq</math> 0.998 or coefficient of determination (r<sup>2</sup>) <math>\geq</math> 0.999</p> <p><b>Specificity</b></p> <p>Detected analysis of interest must not exceed 100% of the mean MCL, or have a <math>\pm</math> SN <math>\geq</math> 10.</p> <p><b>LOD</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>LOQ</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>Recovery value of 80.0% - 120.0% when compared to fresh standard solutions.</p>	
<p><b>Selection of an Appropriate Extraction Solution</b></p> <ul style="list-style-type: none"><li>• Decision should be based on the solubility of the residue</li><li>• Typical extraction solutions utilized include alcohols, waters, buffers or combinations of the three solutions</li></ul>		<p><b>Linearity</b></p> <p>Performed using a minimum of 5 concentration levels, each injected in duplicate. Linearity may range from the limit of quantitation up to 200% of the MCL.</p> <p><b>Specificity</b></p> <p>Typically both swabs and surfaces are evaluated to determine if interferences with the compound of interest are present.</p> <p><b>LOD</b></p> <p>Prepare standard solutions at the estimated LOD (3 preparations) and analyze.</p> <p><b>LOQ</b></p> <p>Prepare standard solutions at the estimated LOQ (3 preparations injected in duplicate) and analyze.</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>The ability of a standard or sample preparation solution to meet method specifications over time.</p>		<p><b>Linearity</b></p> <p>Correlation coefficient (r) <math>\geq</math> 0.998 or coefficient of determination (r<sup>2</sup>) <math>\geq</math> 0.999</p> <p><b>Specificity</b></p> <p>Detected analysis of interest must not exceed 100% of the mean MCL, or have a <math>\pm</math> SN <math>\geq</math> 10.</p> <p><b>LOD</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>LOQ</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>Recovery value of 80.0% - 120.0% when compared to fresh standard solutions.</p>			
<p><b>Rinse or Swab Sampling?</b></p> <p>The sampling technique selected must be capable of quantitatively determining the amount of residual material on the manufacturing equipment. Many considerations go into making the decision on the sampling technique including:</p> <ul style="list-style-type: none"><li>• Ease of access for sampling</li><li>• Size of equipment</li><li>• Solubility characteristics of compound of interest</li></ul> <p><b>Swabbing</b></p> <ul style="list-style-type: none"><li>• Applicable for small surface areas, and difficult to reach areas where traditional swabbing procedures may be difficult</li><li>• Surface should be rinsed long enough to ensure complete coverage and sufficient removal of the target residue</li><li>• More simplistic than swabbing procedures</li><li>• Sample is generally collected from final water rinse of equipment</li></ul> <p><b>Swabbing</b></p> <ul style="list-style-type: none"><li>• Ideal for residues not easily removed with water rinsing</li><li>• Swabs can physically remove insoluble residues</li><li>• Chosen swabs must exhibit the following characteristics:<ul style="list-style-type: none"><li>1) Ability to recover desired residue from given surface</li><li>2) Ability to release the residue into an extraction solution for analysis</li><li>3) Must not contribute excessive interference or background against analysis</li></ul></li></ul>		<p><b>Linearity</b></p> <p>Performed using a minimum of 5 concentration levels, each injected in duplicate. Linearity may range from the limit of quantitation up to 200% of the MCL.</p> <p><b>Specificity</b></p> <p>Typically both swabs and surfaces are evaluated to determine if interferences with the compound of interest are present.</p> <p><b>LOD</b></p> <p>Prepare standard solutions at the estimated LOD (3 preparations) and analyze.</p> <p><b>LOQ</b></p> <p>Prepare standard solutions at the estimated LOQ (3 preparations injected in duplicate) and analyze.</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>The ability of a standard or sample preparation solution to meet method specifications over time.</p>		<p><b>Linearity</b></p> <p>Correlation coefficient (r) <math>\geq</math> 0.998 or coefficient of determination (r<sup>2</sup>) <math>\geq</math> 0.999</p> <p><b>Specificity</b></p> <p>Detected analysis of interest must not exceed 100% of the mean MCL, or have a <math>\pm</math> SN <math>\geq</math> 10.</p> <p><b>LOD</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>LOQ</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>Recovery value of 80.0% - 120.0% when compared to fresh standard solutions.</p>			
<p><b>Swabbing Techniques</b></p> <ul style="list-style-type: none"><li>• Prior to swabbing, swabs are soaked for a few minutes in a vial with the extraction solution to saturate the swab head</li><li>• Excess solution is removed from the swab head by gently pressing the head on the inside of the vial</li><li>• The prepared swabs are rolled with the swab appropriate area using various swabbing patterns</li><li>• Common patterns include partially overlapping parallel swabbing strokes in one direction or back-and-forth. The swab head is then flipped to the other side, and the same pattern is repeated at right angles to the first (see Figure 1.A).</li><li>• Another method is to use a figure-eight swabbing pattern, alternating swabbing strokes in opposite directions, ensuring the swab head never leaves the surface being evaluated. (see Figure 1.B).</li><li>• The swab head is placed back into the vial by dipping the handle above the head with a clean vial</li><li>• One swab may be sufficient to remove residue, but a second or even a third swab can be used to repeat the swabbing pattern to increase the recovery of the residue.</li><li>• Depending on the results used to determine the recovery factor, it may be advantageous to use a dry swab after the wet swab to ensure any remaining solution on the coupon is collected.</li></ul>		<p><b>Linearity</b></p> <p>Performed using a minimum of 5 concentration levels, each injected in duplicate. Linearity may range from the limit of quantitation up to 200% of the MCL.</p> <p><b>Specificity</b></p> <p>Typically both swabs and surfaces are evaluated to determine if interferences with the compound of interest are present.</p> <p><b>LOD</b></p> <p>Prepare standard solutions at the estimated LOD (3 preparations) and analyze.</p> <p><b>LOQ</b></p> <p>Prepare standard solutions at the estimated LOQ (3 preparations injected in duplicate) and analyze.</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>The ability of a standard or sample preparation solution to meet method specifications over time.</p>		<p><b>Linearity</b></p> <p>Correlation coefficient (r) <math>\geq</math> 0.998 or coefficient of determination (r<sup>2</sup>) <math>\geq</math> 0.999</p> <p><b>Specificity</b></p> <p>Detected analysis of interest must not exceed 100% of the mean MCL, or have a <math>\pm</math> SN <math>\geq</math> 10.</p> <p><b>LOD</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>LOQ</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>Recovery value of 80.0% - 120.0% when compared to fresh standard solutions.</p>			
<p><b>Is TOC the Answer?</b></p> <p>Historically, if a compound was rendered by the Merck Index to be only slightly water soluble, everyone immediately ruled out the option of TOC. Lancaster Laboratories has developed methods utilizing TOC for compounds that are not fully water soluble. In order to be TOC compliant, compounds only need to be slightly soluble in water (for example, if a maximum contamination limit (MCL) is set to 10 ppm, the compound only needs to be soluble at 11 ppm in water to be measured in the range of interest). This is a very interesting concept that is growing in popularity across the industry.</p> <p><b>Here are some key factors to evaluate when determining if TOC is the correct approach for a cleaning validation:</b></p> <ul style="list-style-type: none"><li>• Is the carbon content of the sample great enough to be detected once the sample is appropriately diluted?</li><li>• Is the compounds solubility in water above the concentration of the desired MCL?</li><li>• Does the method need to be compound specific?</li><li>• Can acceptable surface recoveries be achieved by rinsing with water or swabbing with water, a weak acid or base?</li></ul>		<p><b>Linearity</b></p> <p>Performed using a minimum of 5 concentration levels, each injected in duplicate. Linearity may range from the limit of quantitation up to 200% of the MCL.</p> <p><b>Specificity</b></p> <p>Typically both swabs and surfaces are evaluated to determine if interferences with the compound of interest are present.</p> <p><b>LOD</b></p> <p>Prepare standard solutions at the estimated LOD (3 preparations) and analyze.</p> <p><b>LOQ</b></p> <p>Prepare standard solutions at the estimated LOQ (3 preparations injected in duplicate) and analyze.</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>The ability of a standard or sample preparation solution to meet method specifications over time.</p>		<p><b>Linearity</b></p> <p>Correlation coefficient (r) <math>\geq</math> 0.998 or coefficient of determination (r<sup>2</sup>) <math>\geq</math> 0.999</p> <p><b>Specificity</b></p> <p>Detected analysis of interest must not exceed 100% of the mean MCL, or have a <math>\pm</math> SN <math>\geq</math> 10.</p> <p><b>LOD</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>LOQ</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>Recovery value of 80.0% - 120.0% when compared to fresh standard solutions.</p>			
<p><b>Application of Correction Factors</b></p> <p>In some cleaning validation studies, it may be determined that not all the residue on the surface can be removed, thus resulting in low recoveries. To increase these results, the application of recovery factors may be necessary. The following is a list of considerations one should evaluate if correction factors are deemed appropriate:</p> <ul style="list-style-type: none"><li>• Recovery factors are usually not applied if recovery results are above 70%, however, there is no standard limit</li><li>• Recovery factors must be set under sound scientific justification</li><li>• Recovery factors should not be used if recoveries are too low (For example, if recoveries are consistently around 10%, a 90% recovery factor would not be appropriate.)</li><li>• Recovery factors need to be set prior to or during validation – not during routine monitoring</li><li>• All results used to determine the recovery factor need to be consistent and reproducible.</li><li>• Method optimization should always be explored as an alternative prior to using recovery factors</li></ul>		<p><b>Linearity</b></p> <p>Performed using a minimum of 5 concentration levels, each injected in duplicate. Linearity may range from the limit of quantitation up to 200% of the MCL.</p> <p><b>Specificity</b></p> <p>Typically both swabs and surfaces are evaluated to determine if interferences with the compound of interest are present.</p> <p><b>LOD</b></p> <p>Prepare standard solutions at the estimated LOD (3 preparations) and analyze.</p> <p><b>LOQ</b></p> <p>Prepare standard solutions at the estimated LOQ (3 preparations injected in duplicate) and analyze.</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>The ability of a standard or sample preparation solution to meet method specifications over time.</p>		<p><b>Linearity</b></p> <p>Correlation coefficient (r) <math>\geq</math> 0.998 or coefficient of determination (r<sup>2</sup>) <math>\geq</math> 0.999</p> <p><b>Specificity</b></p> <p>Detected analysis of interest must not exceed 100% of the mean MCL, or have a <math>\pm</math> SN <math>\geq</math> 10.</p> <p><b>LOD</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>LOQ</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>Recovery value of 80.0% - 120.0% when compared to fresh standard solutions.</p>			
<p><b>How to Handle Failing Data</b></p> <p>The best way to approach this issue is to address it before it is an issue.</p> <ul style="list-style-type: none"><li>• Dedicate a section of the cleaning validation master plan or protocol to handling failing data</li><li>• Lay out a step by step investigation of the results in advance, thus not influencing decisions made on instance by instance circumstances</li><li>• Regulatory agencies like to see that failing results were handled in a consistent and systematic manner</li><li>• All data not meeting acceptance criteria should be handled as a deviation where data should be first verified, resolved, and finally approved</li><li>• Samples may need to be retested or more samples may need to be collected to verify outlying results</li><li>• Modification to the method, protocol, SOP or master plan may need to be entertained if results indicate that a criteria or limit is not attainable</li><li>• All of these scenarios should be investigated during the feasibility/method development/validation stage of the cleaning validation study</li></ul>		<p><b>Linearity</b></p> <p>Performed using a minimum of 5 concentration levels, each injected in duplicate. Linearity may range from the limit of quantitation up to 200% of the MCL.</p> <p><b>Specificity</b></p> <p>Typically both swabs and surfaces are evaluated to determine if interferences with the compound of interest are present.</p> <p><b>LOD</b></p> <p>Prepare standard solutions at the estimated LOD (3 preparations) and analyze.</p> <p><b>LOQ</b></p> <p>Prepare standard solutions at the estimated LOQ (3 preparations injected in duplicate) and analyze.</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>The ability of a standard or sample preparation solution to meet method specifications over time.</p>		<p><b>Linearity</b></p> <p>Correlation coefficient (r) <math>\geq</math> 0.998 or coefficient of determination (r<sup>2</sup>) <math>\geq</math> 0.999</p> <p><b>Specificity</b></p> <p>Detected analysis of interest must not exceed 100% of the mean MCL, or have a <math>\pm</math> SN <math>\geq</math> 10.</p> <p><b>LOD</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>LOQ</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>Recovery value of 80.0% - 120.0% when compared to fresh standard solutions.</p>			
<p><b>When do you need to Revalidate a Cleaning Validation Method?</b></p> <ul style="list-style-type: none"><li>• Cleaning validation procedures should be revalidated when the equipment train of the manufacturing process is changed. Possible changes in the equipment train include the surface type utilized and/or surface area, which can lead to the establishment of a new maximum contamination limit. Usually a full validation can be avoided and only certain elements of the cleaning validation need to be revalidated. If the new limit is within the previously established linear range, only surface recoveries bracketing the new limit and surface residue specificity would need to be revalidated. These same two elements would need to be revalidated if a surface type were changed. If the new limit is outside the previously established linear range, linearity would need to be extended above or below the new limit and swab recovery, surface recovery, and surface residue specificity would need to be revalidated. For a new limit below the established linear range, a new standard concentration at this level may be recommended. However, if the existing method is not linear through the new level, a new standard concentration would be necessary, and this would require full revalidation.</li><li>• Other possible but less likely reasons to revalidate swab recovery, surface recovery, and surface specificity would be a change in the type of swab or swabbing pattern. For a change in swab type, swab specificity would also need to be revalidated. For any of the previously listed changes, elements that would not require revalidation are limits of detection and quantitation and linearity.</li><li>• The prior revalidation discussion assumes that the validated method was for swab samples and not rinse samples. For rinse samples, validation elements involving swabs and surfaces would not need conducted. Additionally, any changes in the synthesis of the drug substance, changes in the composition of the finished product, or changes in the analytical procedure would require revalidation according to ICH guidance.</li></ul>		<p><b>Linearity</b></p> <p>Performed using a minimum of 5 concentration levels, each injected in duplicate. Linearity may range from the limit of quantitation up to 200% of the MCL.</p> <p><b>Specificity</b></p> <p>Typically both swabs and surfaces are evaluated to determine if interferences with the compound of interest are present.</p> <p><b>LOD</b></p> <p>Prepare standard solutions at the estimated LOD (3 preparations) and analyze.</p> <p><b>LOQ</b></p> <p>Prepare standard solutions at the estimated LOQ (3 preparations injected in duplicate) and analyze.</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>The ability of a standard or sample preparation solution to meet method specifications over time.</p>		<p><b>Linearity</b></p> <p>Correlation coefficient (r) <math>\geq</math> 0.998 or coefficient of determination (r<sup>2</sup>) <math>\geq</math> 0.999</p> <p><b>Specificity</b></p> <p>Detected analysis of interest must not exceed 100% of the mean MCL, or have a <math>\pm</math> SN <math>\geq</math> 10.</p> <p><b>LOD</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>LOQ</b></p> <p>Analysis of interest must have a recovery of 75% - 125% of theoretical, with an RSD of <math>\leq</math> 20% (<math>\pm</math> SN of 10).</p> <p><b>Robustness</b></p> <p>The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small, but deliberate variations in method parameters and provides an indication of its reliability during normal usage.</p> <p><b>Stability</b></p> <p>Recovery value of 80.0% - 120.0% when compared to fresh standard solutions.</p>			



## CONCLUSIONS

Establishment of an appropriate cleaning validation platform is critical in any manufacturing process. There are many options to choose from when establishing the cleaning program and care must be taken to ensure that the sampling technique, and analytical monitoring methodology selected, along with establishing the appropriate limit (MCL) and compliance of the approach meet the requirements and the intent of the sample program.

- References:
1. Robeson Bruce Presentation I/V Chicago, (July 2007)
  2. FDA, Guide to Inspections Validation of Cleaning Processes (1993)
  3. D.A. LeBlanc, "Establishing Scientifically Justified Acceptance Criteria for Cleaning Validation of Finished Drug Products," Pharm. Technol. October 1998.
  4. Cleaning Validation and Critical Cleaning Processes, Dohier Research Works
  5. TOC Applications in Pharmaceutical Cleaning Validation, Jon Younkun GE Analytical Instruments, General Electric Powerpoint presentation
  6. Pharmaceutical Technology (April 2009), Using Visible Residue Limits for Introducing New Products into a Pharmaceutical Research Facility
  7. ICH Harmonised Tripartite Guideline, Validation of Analytical Procedures: Text and Methodology, Q2(R1), November 2005

Figure 1.  
There are a variety of patterns that can be followed when using swabs to clean

