Certified Reference Values





Production date: 03-11-2025 Batch: 45 - 2025

		Unhomogenised past. milk	Unhomogenised past. cream										
	Туре	Buttermilk	6,5%	8%	9,5%	11%	12%	14%	18%	32%	38%	46%	50%
Fat (g/100 g) ≤8%+9,5%: ISO 1211/IDF 1:2010 ≥9%: ISO 2450/IDF 16: 2008		0,683	6,692	8,193	9,610	10,807	11,724	13,391	17,595	31,618	38,399	44,778	49,255
Protein (g/ 100 g) ISO 8968-1/IDF 20-1:2014		3,440	3,345	3,302	3,268	3,222	3,549	3,491	3,301	2,636	2,322	2,055	1,849
Lactose, Mono (g/100 g) ISO 5765-2/IDF 79:2002											2,910		
Dry Matter (g/100 g) ISO 6731/IDF 21:2010				16,306	17,600		20,287	21,791	25,596	38,078	44,160	49,982	54,325

** Increase uncertainty (Ux2)

Version 1

Date: 06-11-2025

The results apply only the received batch
The results may not be reproduced, except in their entirety, without written permission from the Laboratory.

Information, if any:

Informations according to the certified reference material:

Description

Bovine milk from Danish Dairies. CRM has preservative added.

The certified values are calulated based on results from chemical analysis.

Because of the short validy period, reference analyses is carried out at one laboratory only, but in several determinations and batches for statistical robustness.

Homogenity iscontrolled on data from both chemical and FTIR analysis.

Stability has been calculated for a representative batch.

*Following changes to the above methods have been made.

Fat: The drying is with a constant time instead of constant weight.

Protein: The destruction does not take place at lower temperatures, but starts at 420°C.

Dry matter: Single-use vials without lids are used instead of vials with lid, and the samples are not pre-dried before placed in the oven.

The Samples are calculated with following expanded uncertanity U%

	U%
Fat	1,5%
Protein	1,5%
Lactose	2,0%
Dry Matter	1,0%

Instruction for handling and use:

Place the CRM at room temperature for 60 min, then heat in water bath to 38°C for 20 min.

Gently invert the CRM 20 times and measure immediately on FTIR equipment

Samples can't be reused. Once heated, please discard any remains after measurements have finished.

Sample size: Approx. 50 mL.

 $\frac{Period\ of\ validity:}{14\ days\ at\ 4^{\circ}C\ \pm\ 2^{\circ}C\ from\ production\ date}.$

Storage information:

Store at 4°C ± 2°C until use.

°C.