

ABN 50 005 085 521

SAMPLE CONTAINER ORDER FORM

Client Name:	Date Requested:
Site Name/Job #	
Delivery Address:	Date Required:
Contact Name:	Sent by:
Phone Number :	Checked by:
Eurofins mgt contact:	
Service required	OVERNIGHT FIRST CLASS* LOCAL COURIER
Please circle:	COLLECT SAME DAY* Eurofins mgt COURIER
	INTERNATIONAL FREIGHT COST TO BE ON-CHARGED? Y / N

DELIVERY: LOCAL (Auckland/Wellington/Christchurch) allow two days. REGIONAL allow five days. *For urgent deliveries contact Client Services.

Colour	Qty Req	Description	Typical Tests	No/sample
Green		Single jars – 250mL	All soil analysis	1
Green		Boxes of 48 jars – 250mL	All soil analysis	1
Grey		Plastic Bag – press sealable	Acid sulphate soils	1
Dark Blue		Plastic Bag – press sealable	Asbestos	1
Green		1000mL Plastic	pH, BOD, SS, TDS, EC, Alkalinity, NO3	1
Orange		1000mL Glass	SVOC (TRH, PAH, PCDD/F, PBDEs), chlorophyll-a	1#
Orange		200mL Glass	Semi-volatiles (TRH, PAH etc), Glyphosate	1#
Purple		250mL Glass – Preserved with H ₂ SO ₄	Oil & Grease	2
Purple		125mL Glass – Preserved with H ₂ SO ₄	TOC, COD	1
Pink		40mL VOA Vials – Preserved with Sodium Bisulfate	BTEX, MAH, VOC, Volatiles, Gases (methane etc)	2
Black / Grey		4L Amber - Preserved with CuSO ₄	Dust Gauges	1
Purple		250mL Plastic – Nutrients – Preserved with H ₂ SO ₄	Ammonia, Total phenolics, PO ₄ , TKN	1
Black / Grey		200mL Plastic – No Preservative	Metals, dissolved - not field filtered (not recommended)	1
Red		200mL Plastic – Metals – Preserved with Nitric acid	Metals, dissolved - field filtered or metals, total	1
Mauve		200mL Plastic – Ferrous – Preserved with HCl	Ferrous Iron	1
Blue		125mL Plastic - Preserved with NaOH	Cyanide, Chromium (hexavalent)	1
Yellow		125mL Plastic - Preserved with Zinc Acetate + NaOH	Sulphide	1
Grey		125mL Plastic - Preserved with EDTA	Sulphite	1
Grey		Sterile 500mL plastic - Preserved with Na ₂ S ₂ O ₃	E. Coli, Plate counts, Coliforms	1
Grey		Sterile 500mL plastic - Preserved with Na ₂ S ₂ O ₃	Salmonella, Campylobacter	2
Grey		Sterile 120mL plastic - Preserved with Na ₂ S ₂ O ₃	Legionella	1
Grey		Sterile 120mL plastic soil jar	Microbiological analysis in soil	1
Dark Blue		500mL PET	PFOS/PFOA/AFFF in water	1
Dark Blue		100mL white polypropylene jar	PFOS/PFOA/AFFF in soil	1

Other	Quantity	Packing	Quantity
		Large Esky	
		Medium Esky	
		Small Esky (Plastic)	
		Small Esky (Foam) for Soil or Water	
		Ice bricks	
		Bubble wrap bags - for glass bottles	

Where low level detection required refer to your Analytical Service Manager regarding sample volumes

NB: We recommend that you provide additional sample on the 1st, 11th, 21st, etc sample for performance of Duplicates / Matrix Spikes. (Note however that Matrix spike determinations are not appropriate for all tests).