



Instrumentation - Barcelona, Spain

Liquid Chromatography

Waters HPLC	1
Waters HUPLC	3
Waters UPLC	1
Agilent HPLC	3
Thermo Dionex HPLC	3
HPLC Detectors	
Waters Ultraviolet (UV)	2
Waters Photodiode Array (PDA)	3
Agilent Photodiode Array (PDA)	3
Waters Refractive Index (RI)	1
Agilent Refractive Index (RI)	1
Agilent Evaporative Light Scattering Detector (ELSD)	1
Waters Fluorescence (FL)	1
Dionex Conductivity Detector (Ion Chromatography)	2
Dionex Electrochemical Detector (Ion Chromatography)	3
Thermo Dionex Ultraviolet	2
Dionex Pulsed Amperometric Detector	3
Post Column Derivation Module	2
Waters GPC Empower Software	1

Gas Chromatography

Agilent/Bruker/Konik GC	8
GC Detectors	
Flame Ionisation Detector (FID)	9
Electron Capture Detector (ECD)	2
Thermal Conductivity Detector (TCD)	3
Flame Photometric Detector (PFPD)	1
Nitrogen-Phosphorus Detector (NPD)	2
Liquid Autosamplers	6
Headspace Autosamplers	2

Mass Spectrometry

Agilent GC/MS Triple Quad 7000A with Liquid Autosampler	1
Agilent GC/MS Triple Quad 8890 with Liquid Autosampler	1
Agilent GC/MS Triple Quad 8890 with Combipal Headspace and SPME	1
Thermo Fischer LC-MS Orbitrap Exploris 120	1
Thermo Fischer LC-MS Orbitrap Exploris 240	1
Waters LC-MS Xevo TQ XS	1

Elemental Analysis

Perkin Elmer Graphite Furnace Atomic Absorption Spectrophotometer	2
Perkin Elmer Flame (Acetylene) Atomic Absorption Spectrophotometer	1
Perkin Elmer Flame (Acetylene-Nitrous Oxide) Atomic Absorption Spectrophotometer	1
Perkin Elmer Hydride Generation Atomic Absorption Spectrophotometer	1



Perkin Elmer Inductively Coupled Plasma (ICP-OES) Spectrophotometers	2
Perkin Elmer Ultrasonic Nebuliser ICP Autosamplers	1
Perkin Elmer Cold Vapour and Hydride (Hg) Atomic Absorption Spectrophotometer	1
Perkin Elmer ICP-MS	1
Milestone Ethos Microwaves	2
Milestone Ultrawave ECR	1

Dissolution

PharmaTest Dissolution Bath	4
PharmaTest Disaggregation Bath	1

Spectrophotometry

Perkin Elmer Fluorescence Spectrophotometer	1
Jasco IR Spectrophotometer	2
Aligent UV/Vis Spectrophotometer Cary 60	2

Biochemistry

Bio Tek Micro Plate Reader (ELISA)	1
------------------------------------	---

Cell Culture

Microtome	1
Tissue Embedding System	1

Microbiology

Microbiology Laboratory with Controlled Access (Clean Room Class D)	1
Sterilisation and Depyrogenation Devices (Autoclaves and Oven)	3
Biohazard Hoods	2
Microbiology Cabinets	6
Millipore Filtration Systems	6
Geringe La Calhene Isolator 4 Gloves (Sterility Test)	1
Steritest Equinox Pump (Isolator Filtration System)	2
Active Air Sampler	1
SKAN Pure 2 Isolator 4 gloves (Sterility Test)	1
Steritest Symbio Pump (Isolator Filtration System)	
Active Air Sampler	
Incubators	14
CO ₂ Incubators	2
Refrigerators/Freezers/UltraFreezers	14
Water Baths	8
Ultrasonic Baths	4
Microscopes	2
Colony Counter Digitals S	6
Biotech Synergy Mx with Spectrophotometer and Luminometer Module	1
Biotech/Agilent Kinetic System for Performing Kinetic Chromogenic and turbidimetric Endotoxin Tests	2
Novasina Water Activity	2
Aqualab Water Activity	1
Miele Washer Disinfectors	1



Particle Counter PAMAS SBSS and SVSS (viscous samples)	2
Gas analyser Dräger Multitest	1
Gas sampler Pinnocchio Super 2	1

Stability

CTS/Dycometal Climatic Chambers	8
25°C±2°C/60%HR ±5%HR	40 m ³
25°C±2°C/40%HR ±5%HR	2,6 m ³
30°C±2°C/65%HR ±5%HR	8 m ³
30°C±2°C/75%HR ±5%HR	1,1 m ³
30°C±2°C/75%HR ±5%HR	1,1 m ³ + 1,2 m ³
40°C±2°C/75%HR ±5%HR	7,5 m ³
5°C±3°C	11,3 m ³
Custom Conditions	1,3 m ³

Thermal Analysis

Buchi Melting and Boiling Point Analyser	1
--	---

Facility Monitoring (Water, Gas, Clean Rooms)

Total Organic Carbon Analyser (TOC)	2
Air Samplers	5
Conductivity Meters	2
pH Meters	8
Portable Particle Counter (0.3 – 5 µm)	1

Miscellaneous

Viscosimeters	
Brookfield Viscosimeter LU/ULA	1
Brookfield Viscosimeter HB	1
Brookfield Viscosimeter RVDV-II	1
Brookfield Recirculation Bath	1
Cole-Palmer Masterflex® Peristaltic Pumps	2
Peristaltic Pumps	3
Rudolph Polarimeter: 589nm-546nm-633nm-426nm-405nm-365nm	1
Fiske Osmometer	1
Advanced Technologies Osmometer Osmotech	1
HIAC/Beckman Liquid Particle Counter	1
PharmaTest Durometer	1
Bausch & Lomb Refractometer	1
Methrom Volumetric Automatic System	1
Methrom Volumetric Karl Fischer (KF) System	2
Methrom Coulometric Karl Fischer (KF) System	2
Methrom Electronic Titrations System	1
Vaccum Oven	1
Vaccum Oven (up to 7mBar)	1
Muffle Furnace	2
Fluorescence Visualiser for Thin-Layer Chromatography (TLC)	1



Microbalance	1
Analytical Balances	8
Scaltec Moisture Balance (IR Balance)	1
Friabilimeter	1
Sieving Machine	1
Buchi Distillation Unit K-355 (Kjedahl and Devarda Method)	1
Buchi Sample Speed Digester	1
Hazardous Material Cabinet	5
Turbovap Evaporator	1

Metrology

MySirius (Facility Monitoring System - Temperature and Humidity)	1
Full-Time Metrologist On Site	