



Instrumentation – Sweden

Liquid Chromatography

Shimadzu HPLC LC10	1
Shimadzu HPLC LC20	4
Shimadzu HPLC LC30	2
HPLC Detectors	
Ultraviolet (UV)	4
Photodiode Array (PDA)	3
Refractive Index (RI)	3
Evaporative Light Scattering (ELSD)	2
Fluorescence (FL)	1
Conductivity (CCD)	1

Gas Chromatography

Shimadzu GC2010	2
Shimadzu GC2010+Headspace	2
Shimadzu GC2014 (FID, TCD)	1

Mass Spectrometry

Shimadzu LCMS 2020 (Single Q)	1
Shimadzu QP2010 Ultra GCMS+Headspace	1

Dissolution

Sotax Dissolution AT7	2
Sotax Disintegration DT2	1
Sotax HT1 Hardness Tester	1
Sotax F1 Friabilator	1
Pharma Test Suppository Disintegration Tester	1

Spectrophotometry

Shimadzu UV 2401	1
Shimadzu IRAffinity +ATR	1
Jasco Polarimeter P-2000	1

Stability

Termaks 400L	7
Nüve 620	5

Miscellaneous

Metrohm Micro KF 832	1
Metrohm KF 841	1
Shimadzu TOC-V	2
Anton Paar Density Meter	1
HACH Turbidimeter 2100P	1
Mettler Refractometer	1
MYR Viscosity Meter	2
SRS MPA 161 Melting Point Analyzer	1
WTW Conductivity Meter	2
Metrohm 905 Titrand	1