



Instrumentation – Livingston, UK

Liquid Chromatography

HPLC Detectors

Agilent 1100 UV Detectors	4
Agilent 1200 UV Detectors	2
Agilent 1260 UV Detectors	7
Agilent 1100 Diode Array Detectors	4
Agilent 1100 RI Detector	1
Agilent 1260 Diode Array Detectors	3
Agilent 1260 RI Detectors	2
Dionex Ultimate 3000, Diode Array Detector	1

Gas Chromatography

Hewlett Packard, HP6890 Series/G1530A, GC-FID with Liquid Autosampler	1
Hewlett Packard, HP6890 Series/G1530A, GC-FID and GC-TCD with Liquid Autosampler	1
Agilent 6850 GC-FID with Liquid Autosampler	1
Agilent 7890 GC-FID with Liquid and HS Autosamplers	2
Agilent 7890 GC-FID with Liquid Autosampler	3

Mass Spectrometry

Agilent 7890 GC-MS with Liquid Autosampler	1
--	---

Elemental Analysis

Thermo AA Spectrophotometer	1
Agilent 7700 ICP-MS	1

Dissolution

Caleva Dissolution Apparatus 10ST+ (Type 1, 2)	4
Copley Dissolution Apparatus DIS6000	2
Sotax Dissolution Apparatus AT7 (Type 2)	4

Spectrophotometry

Kontron Fluorescence Spectrophotometer	1
Perkin Elmer FT-IR Spectrophotometer	1
Shimadzu UV-1800 UV-Vis Spectrophotometer	1

Microbiology

Autoclaves	2
Air Samplers	2
Microscopes	2
Clean Room	1
Colony Counter	1
Isolators	2
LAF Cabinets	4



Miscellaneous

Mettler Toledo Automatic Titrator	1
Hanna Conductivity Meters	2
Anton Paar Density Meter	1
Copley Scientific Disintegration Tester	1
Copley Scientific Friability Tester	1
Hardness Testers	2
HACH 9703+ Liquid Particle Counter	1
Mettler Toledo Melting Point Apparatus	1
Microscope (Polarising)	1
Milestone Microwave Digester with TraceCLEAN	1
KF Mettler Toledo (Coulometric)	1
KF Mettler Toledo (Volumetric)	1
Carbolite Muffle Furnaces	2
Vitech Osmometer	1
pH Meters	4
Rudolph Polarimeter	1
Bellingham Refractometer	1
Brookfield Viscometer	1
Bellingham and Stanley Polarimeter	1

Storage

Refrigerators	6
Liebherr Freezer	1
Incubators	4
Walk-in Cold Storage (3mx3m)	2

Metrology

Anemometers	2
Mitutoyo Ltd Digimatic Vernier Caliper	1
Digital Lux Meter	1
TPI Digital Manometer	1
RS Components Digital Photo Tachometers	2
Digi-Pas Digital Protractor	1
Kinesis HPLC Liquid Flow Meters	2

Instrumentation – Camberley, UK

Stability Storage

Walk-in Stability Room (25°C/60%RH)	1
Stability Cabinets (2°C up to 50°C/75%RH)	9
Reach-in Stability Rooms (2°C up to 40°C/75%RH)	11
Stability Freezers (-80°C and -20°C)	4
Photostability Cabinets	2
Stability Ovens (50°C to 60°C)	3

Miscellaneous

Eurotherm Temperature and Humidity Monitoring Systems	2
Eurotherm Portable Temperature and Humidity Monitoring System	1
Analytical Balances	3
Balances	3