



Accredited Laboratory

A2LA has accredited

MET LABORATORIES, INC.

Austin, TX

for technical competence in the field of

Electrical Testing

This laboratory is accredited in accordance with the recognized International Standard ISO/IEC 17025:2005 *General requirements for the competence of testing and calibration laboratories*. This accreditation demonstrates technical competence for a defined scope and the operation of a laboratory quality management system (refer to joint ISO-ILAC-IAF Communiqué dated 8 January 2009).



Presented this 27th day of June 2017.

A handwritten signature in black ink, appearing to read 'L. Sen', written over a horizontal line.

President and CEO
For the Accreditation Council
Certificate Number 0591.06
Valid to March 31, 2019

For the tests to which this accreditation applies, please refer to the laboratory's Electrical Scope of Accreditation.