



Internet of Things

Accredited test laboratories, certification bodies and competence centers provide services for placing products on global markets

- Wireless alliance services (Bluetooth SIG, Qi wireless charging, PTCRB, etc.)
- Regulatory compliance test and certification (CE-marking, automotive E-marking, US FCC, etc.)
- Manufacturer compliance tests (MBN, VW TL, JLR, SAE, etc.)
- International homologation for radio technologies (Chinese SRRC & NAL, Japanese MIC, etc.)
- International homologation for safety and EMC (Chinese CCC, Japanese PSE, Korean KCC, etc.)
- Reliability tests environment simulation (mechanical, climatic, etc.)
- Interoperability & functional testing (MirrorLink, Serial ATA, USB, DLNA, etc.)
- Test automation (Smartphone app, etc.)



The competence of more than 15 test laboratories in one partner

- | | | |
|---------------------------|----------------------------|---------------------------|
| • Belgium (Diepenbeek) | • Italy (Torino) | • Switzerland (Bern) |
| • Germany (Reichenwalde) | • Scotland (Grangemouth) | • Switzerland (Rossens) |
| • Great Britain (York) | • Switzerland (Fehraltorf) | • Switzerland (Rorschach) |
| • Great Britain (Bristol) | • Switzerland (Zürich) | • China (Shanghai) |
| • Great Britain (Leeds) | • Switzerland (Neuenhof) | • Thailand (Bangkok) |

Contact our expert team: **IoT@eurofins.com**