



Instrumentation – Sydney, Australia

Gas Chromatography

Agilent GC	1
Agilent Headspace Sampler	1

Spectrophotometry

Shimadzu UV 160	1
-----------------	---

Molecular Biology/Cell Biology/Virology

Hitachi 3500 Genetic Analysers Model ABI3500	2
Invitrogen E-Gel i Base E-Gel Safe Imager Real Time Transilluminators	2
ThermoFisher Scientific Centrifuge Heareus Megafuge 8	1
Forma CO ₂ Incubator	1
IGN Flow Co ₂ Incubator	1
Veriti Thermacyclers	2

Microbiology

Biomeriux Vitek/Vitek MS	1
Olympus Microscope	1
Biological Safety Cabinets Class 2	6
Laminar Flow Units	6
ThermoFisher Cytotoxic Drug Safety Cabinets	2
AI Stomacher	1
Millipore Manifolds	2
Water Baths	5
Walk-In Incubators	5
Stand-Alone Incubators	6
Biotex Ultra Micro Plate Readers	2
Haraeus Centrifuge Multifuge	1
Haraeus Centrifuge Fresco	1
Particle Measuring System Particle Counter APSS-2000	1
Millipore Manifold Base	1
Millipore Universal Pump	1
Syngene G-Box HR	1
Eppendorf PCR Thermocycler Mastercycler	1
Labnet Digital Dry Bath Accublock	1
Applied Biosystems Real-Time PCR System StepOne	1
Atherton Autoclave	1
Labec Autoclave	1
Shimadzu Spectrophotometer	1
Lotus Microscope	1
Colony Counters	2

Metrology

Sartorius Analytical Balance	1
Top Pan Balances	3

Stability

Walk-in Cold Room	1
-------------------	---

Miscellaneous

Hanna pH Meter	1
Hanna Dissolved Oxygen Meter	1
Hanna Turbidity and Chlorine Meters	2
Hanna Thermo-Hygrometer	1



Miscellaneous (contd)

Labchem pH Meter	1
Fluorescence Analysis Box	1
Unisonics Ultrasonic Cleaner	1
CCL Biological Safety Cabinet Class 2	1
Novasina Water Activity Meter	1

Instrumentation – Melbourne, Australia

Liquid Chromatography

Agilent 1100 HPLC's	7
Agilent 1200 HPLC's	10
Agilent 1260 HPLC's	2
Agilent 1290 HPLC's	2
Metrohm, 960 Professional IC Vario Ion Chromatography System	1

Gas Chromatography

Hewlett Packard, HP6890 GC	5
Agilent, 7890B GC	2
Agilent, 8890 GC	1
Agilent 7697A Headspace Samplers	2

Mass Spectrometry

Agilent 5975 MSD (GC)	1
Agilent 6410 Triple Quad MS (LC)	2

Elemental Analysis

Varian 730 ES, ICP-OES	1
Varian Spectra AA-200 Atomic Absorbance Spectrometer	1

Dissolution

Vankel VK7000 Dissolution Baths	4
---------------------------------	---

Spectrophotometry

Mettler Toledo UV7	1
--------------------	---

Microbiology

Associates of Cape Cod, Pyros Kinetix Flex System, PKF64	1
--	---

Stability

5C, 25C/60%RH, 30C/65%RH, 30C/75%RH & 40C/75%RH	5
Photostability Chamber	1

Metrology

Mettler Toledo Analytical Balances	6
Top Pan Balances	3