

STERILE MANUFACTURING

Primary packaging & finished dosage forms

EXPERTISE

- **10 years of experience in Fill & Finish**


Manufacturing process under isolators
Manual and semi-automatic filling lines

- **Authorizations for all DS except sexual hormones, cephalosporin/β-Lactam**

Biotech product, nanocarrier, highly potent product, controlled substance

- **Processing of various small molecules to biologics**

Amino acid, peptide, oligonucleotide, plasmid vector, cytokines, protein, immunoglobulin, interferon, mAb, ADC, inactivated bacterias...



560 000
MFT units filled
without
any positive
results

1000+
batches

SERVICES

Fill & Finish

- Compounding
- Aseptic filling
- Pilot, clinical and small scale commercial batches
- Viscous products

QP Confirmation & Release



STERILE MANUFACTURING

Primary packaging & finished dosage forms



FILL & FINISH CAPACITIES



VIALS

PREFILLED SYRINGES

Small size of bulk volume

Type of packaging	Volume/Units per batch
Bulk	From < 0,5 L & up to 40 L
Vials	Up to 5,000 units (e.g. 2R vials)
Syringes	Up to 5,000 units (e.g. syringes 1 mL long)


Small filling volumes
from 0.5 mL
& up to 50 mL

Other packaging upon request

FACILITIES

FILL & FINISH

- 2,200 m² / 24,000 ft² dedicated to manufacturing
- 480 m² / 5,167 ft² classified C clean rooms
- Storage capacities at RT, 2-8°C, -20°C, -80°C
- Experienced staff in sterile manufacturing
- State-of-the-Art equipment per industrial standard but designed for low scale batches
- Total waste management strategy (disposable materials & solvents)



ANSM ^{GMP}

PMDA

FDA



MAY 2022

