



Instrumentation – Denmark

Liquid Chromatography

Agilent 1100 HPLC	1
Waters Alliance HPLC	1
Dionex Ultimate HPLC	2
HPLC Detectors	
Agilent 1100 Variable Wavelength Detectors (UV)	2
Agilent 1100 Refractive Index Detector G1362A (RI)	1
Dionex Refractive Index Detector Shodes RI-101 (RI)	1
Dionex Photodiode Array Detector UVD340U (PDA)	1
Dionex VWD-3100 Detector	1
Dionex Diode Array Detector DAD-3000 (DAD)	2

Gas Chromatography

Agilent 6890A GC-FID	2
Perkin Elmer Headspace Sampler TurboMatrix 40 HS	2

Ion Chromatography

Dionex IC ICS-1500	1
Dionex IC ICS-1600	1

Metal Analysis

Varian GF-AAS AA240Z GTA 120	1
Varian FL-AAS AA240FS	1
Thermo ICP-AES ICAP 6300 and 7400	2
Milestone Microwave Oven Ethos	1
Milestone Ultrawave	1
Anton Paar Microwave Oven Multitwave	2
Agilent ICP MS 7700x	1
Agilent ICP-MS 7900	3

Dissolution

Distek Dissolution Evolution 6300	1
-----------------------------------	---

Spectrophotometry

Thermo Fischer Spectrofotometer Evolution 60	1
Thermo FT-IR Spectrofotometer Nicolet 380	1

Microbiology

Termarks Incubators	18
Laf Benches	6
Getinge Autoclave	1
Biotek Instruments Microplate Reader	1
Cold Storage Rooms (5° C)	2
AB 3500 XL Genetic Analyzer (MicroSeq)	1
Veriti Thermal Cycler	1


Microbiology (cont.)

Binder Incubators	2
Lab refrigerators	8
Lab Freezer, -80° C	2
Lab Freezer, -20° C	1
Labguard for temperature and humidity monitoring	1

Stability

Sanyo PSC062.XHX.C	3
Binder KBF 720, 25°C/60%RH	2

Miscellaneous

Methrohm pH-meter 713	1
Shimadzu TOC-V CPH	1
Shimadzu TOC-V CSH	1
Radiometer Conductometer CDM83	1
Methrohm Autotitrator 702 SM Titrino	1
Anton Paar Densitometer DMA38	1
B&S Polarimeter ADP220	1
Erweka Disintegrator ZT 3-1 & Type E	1
Mettler Toledo Karl Fisher Titrator DL38	1
Thermo/HAAKE Rotationviscosimeter VT550	1

Metrology

Mettler Toledo AT 261 Delta range	2
Mettler XP205	1
Mettler XP603S DR	1

For equipment details, see VMP.