

Analysis type	Instrument	Number of instruments
Liquid Chromatography	Agilent 1100 HPLC	1
	Waters Alliance HPLC	1
	Dionex Ultimate HPLC	2
	HPLC Detectors	
	Agilent 1100 Variable Wavelength Detectors (UV)	1
	Dionex Refractive Index Detector Shodes RI-101 (RI)	1
	Dionex Photodiode Array Detector UVD340U (PDA)	1
	Dionex VWD-3100 Detector	1
	Dionex Diode Array Detector DAD-3000 (DAD)	1
Gas Chromatography	Agilent 6890A GC-FID	2
	Perkin Elmer Headspace Sampler TurboMatrix 40 HS	2
Ion Chromatography	Dionex IC ICS-1500	1
	Dionex IC ICS-1600	1
Metal Analysis	Thermo ICP-AES ICAP 6300	1
	Thermo ICP-AES ICAP 7400	1
	Milestone Microwave Oven Ethos	1
	Anton Paar Microwave Oven Multitwave	3
	Agilent ICP MS 7700x	1
	Agilent ICP-MS 7900	2
Spectrophotometry	Thermo Fischer Spectrofotometer Evolution 60	1
	Thermo FT-IR Spectrofotometer Nicolet 380	1
Microbiology	Termaks Incubator	18
	Binder Incubator	2
	LAF Benches	6
	Getinge Autoclave	1
	Biotek Instruments Microplate Reader	1
	Cold Storage Rooms (5o C)	2
	AB 3500 XL Genetic Analyzer (MicroSeq)	1
	Veriti Thermal Cycler	2
	TOMY Autoclave	1
	Lab refrigerators	7
	Lab Feezer, -80 C	1
	Lab Freezer, -20 C	2
	Labguard for temperature and humidity monitoring	1

Analysis type	Instrument	Number of instruments
Miscellaneous	Mettler Toledo Seven Compact S220 pH meter	1
	Shimadzu TOC-V CPH	1
	Shimadzu TOC-V CSH	1
	Mettler Toledo Seven Multi Conductometer	1
	Methrohm Autotitrator 702 SM Titrino	1
	Anton Paar Densitometer DMA38	1
	B&S Polarimeter ADP220	1
	Mettler Toledo Karl Fisher Titrator DL18	1
Metrology	Mettler Toledo AT 261 Delta range	1
	Mettler XP205 DR	3
	Mettler XP603S DR	1

Instrument details are in instrument files and eLIMS-BPT Metrology Module