

Type	No.	Fat %	95%Conf.	Protein %	95%Conf.	Dry Matter %	95%Conf.	Lactose	95%Conf.	Casein	95%Conf.
30	R1.5	1.5 (1.2-1.8)	±0.05	3-4	±0.03	10-11.5	±0.06	4.5-5.5	±0.08	-	-
C	R3.0	3.0 (2.5-3.5)	±0.02	3-4	±0.03	11-13	±0.07	4.5-5.5	±0.08	-	-
B	R4.0	4.0 (3.6-4.4)	±0.03	3-4	±0.03	12-14	±0.07	4-5	±0.08	2.6	±0.08
A0	R5.5	5.5 (4.5-6.5)	±0.04	4-5	±0.04	14.5-16.5	±0.09	-	-	-	-
A1	R6.1	6.1 (5.1-7.1)	±0.04	3.5-4.5	±0.04	14-17	±0.09	4-5	±0.08	-	-
6.5%	RU6.5	6.5 (6.0-7.0)	±0.04	3-4	±0.03	-	-	-	-	-	-
8%	RU8.0	8.0 (7.5-8.5)	±0.05	3-4	±0.03	15.5-17	±0.09	-	-	-	-
9.5%	RU9.5	9.5 (9.0-10.0)	±0.06	2.5-4	±0.03	-	-	-	-	-	-
11%	RU11	11 (10.3-11.7)	±0.07	2.5-4	±0.03	-	-	-	-	-	-
Skm	PU0.07	0.07 (0.04-0.10)	±0.02	3-4	±0.03	9-10	±0.05	4.5-5.5	±0.09	-	-
Sweet But	PU0.7	0.7 (0.5-0.9)	±0.04	2.5-4	±0.03	-	-	-	-	-	-
1%	OU1	1 (0.5-1.5)	±0.06	3-4	±0.03	9-11	±0.06	4.5-5.5	±0.08	2.8	±0.08
3%	OU3	3 (2.5-3.5)	±0.02	3-4	±0.03	11-13	±0.07	4.5-5	±0.09	2.7	±0.08
5%	OU5	5 (4.5-5.5)	±0.03	3-4	±0.04	13-14.5	±0.08	4-5	±0.09	2.9	±0.09
7.5%	PC7.5	7.5 (7.0-8.0)	±0.05	2.9-3.6	±0.03	15-17	±0.09	-	-	-	-
9%	PC9	9 (8.5-9.5)	±0.06	2.9-3.6	±0.03	16-18	±0.10	-	-	-	-
12%	PC12	12 (11.3-12.7)	±0.07	2.8-3.4	±0.04	19-21	±0.11	-	-	-	-
14%	PC14	14 (13.4-14.6)	±0.08	2.7-3.4	±0.04	21-22.5	±0.12	-	-	-	-
18%	PC18	18 (17.4-18.6)	±0.10	2.5-3.2	±0.04	24.5-26	±0.14	-	-	-	-
32%	PC32	32 (31.5-32.5)	±0.20	2.0-2.7	±0.03	37-39	±0.21	-	-	-	-
38%	PC38	38 (37.2-38.8)	±0.20	1.9-2.4	±0.03	43-45	±0.24	2.5-3.5	±0.06	-	-
42%	PC42	42 (40.7-43.3)	±0.30	1.7-2.5	±0.03	46-48	±0.27	-	-	-	-
48%	PC48	48.5 (47.5-49.5)	±0.30	1.5-2.0	±0.03	52-54.5	±0.31	-	-	-	-
Milk 0.5	PH0.4	0.4 (0.2-0.6)	±0.03	3-4	±0.03	9-10.5	±0.06	-	-	-	-
Milk 1.5	PH1.5	1.5 (1.3-1.7)	±0.09	3-4	±0.03	10-11.5	±0.06	-	-	-	-
Milk 3.0	PH3.0	3.0 (2.7-3.3)	±0.02	3-4	±0.03	11-12.5	±0.07	-	-	-	-
Milk 3.5	PH3.5	3.5 (3.2-3.8)	±0.02	3-4	±0.03	12-13	±0.07	-	-	-	-
Milk 4.0	PH4.0	4.0 (3.7-4.3)	±0.03	3-4	±0.03	-	-	-	-	-	-
Whey 1	WP0.07	0.07 (0.04-0.10)	±0.03	0.6-1.1	±0.01	4-6	±0.03	3-4	±0.07	-	-
Whey 2	WP0.15	0.15 (0.09-0.21)	±0.03	1.5-3	±0.03	6-8	±0.04	3.5-5	±0.08	-	-
Whey 3	WP0.22	0.22 (0.15-0.29)	±0.05	3-4	±0.03	8.5-10.5	±0.05	4-5.5	±0.08	-	-
Junket	CM3.5	3.5 (3.0-4.0)	±0.03	3-4	±0.03	-	-	-	-	-	-
Butter Milk	CM0.5	0.5 (0.25-0.75)	±0.10	3-4	±0.03	-	-	-	-	-	-