



## *Accredited Laboratory*

A2LA has accredited

### **EUROFINS PRODUCT SERVICE GMBH**

*Reichenwalde, Germany*

for technical competence in the field of

### **Electrical Testing**

This laboratory is accredited in accordance with the recognized International Standard ISO/IEC 17025:2005 *General requirements for the competence of testing and calibration laboratories*. This accreditation demonstrates technical competence for a defined scope and the operation of a laboratory quality management system (refer to joint ISO-ILAC-IAF Communiqué dated 8 January 2009).



Presented this 23<sup>rd</sup> day of August 2016.

A handwritten signature in black ink, written over a horizontal line.

Senior Director of Quality and Communications  
For the Accreditation Council  
Certificate Number 1983.01  
Valid to October 31, 2018  
Revised September 27, 2018

*For the tests to which this accreditation applies, please refer to the laboratory's Electrical Scope of Accreditation.*