

ACCREDITATION CERTIFICATE

No. LA.01.142

National accreditation body hereby certifies that

complies with the requirements of

„Eurofins Labtarna Lietuva“, UAB

LST EN ISO/IEC 17025:2018

legal entity code: 123647492

and is competent to perform this accredited activity:

microbiological tests of water, food, feeding stuffs and environmental samples

The scope of accreditation below is an integral part of this certificate. Locations of the conformity assessment body are specified in the scope of accreditation

Initial accreditation date: **2015-06-15**

Certificate issued / valid since: **2025-06-05**

Version of: **2026-06-12**

Expiry date: **2030-06-04**

Director of the Accreditation Department



TADAS JUODELIS

The certificate may be changed, its validity suspended or withdrawn by the decision of the National accreditation body. Information on the actual data of accreditation certificates may be verified at lasa.lrv.lt



„Eurofins Labtarna Lietuva“, UAB, accredited in accordance with LST EN ISO/IEC 17025:2018

Location of the conformity assessment body:

Kauno st. 1A, LT-01314 Vilnius

Materials or products tested	Component, parameter or characteristic to be tested	Reference number of the document specifying test methods, clause (if relevant)	Techniques, methods and/or equipment used (where appropriate)
Microbiological tests of food and feeding stuffs			
Food, feeding stuffs and environmental samples	Enumeration of microorganisms, 30 °C	LST EN ISO 4833-1:2013 LST EN ISO 4833-1:2013/A1:2022	Counting method. Pour plate technique
	Detection of <i>Salmonella</i> spp.	LST EN ISO 6579-1:2017 except LST EN ISO 6579-1:2017/A1:2020	Detection method. Principle of enrichment and surface inoculation
	Detection of <i>Listeria monocytogenes</i>	LST EN ISO 11290-1:2017	Detection method. Principle of enrichment and surface inoculation
	Enumeration of <i>Enterobacteriaceae</i> , 37 °C	LST EN ISO 21528-2:2017	Counting method. Pour plate technique
	Detection of <i>Salmonella</i> spp.	SVP 5.4 M-50, 2012-02	Polymerase chain reaction (PCR). BAX system
Food stuffs and environmental samples	Enumeration of Mesophilic Lactic acid bacteria	LST ISO 15214:2009	Counting method. Pour plate technique
	Detection of <i>Escherichia coli</i>	LST ISO 7251:2005 LST ISO 7251:2005/A1:2024	Detection method. Principle of inoculation into a liquid medium
	Detection of coliforms	LST ISO 4831:2006	Detection method. Principle of inoculation into a liquid medium
Food and feeding stuffs	Enumeration of coliforms	LST ISO 4832:2006	Counting method. Pour plate technique
	Enumeration of β -glucuronidase-positive <i>Escherichia coli</i>	LST ISO 16649-2:2002	Counting method. Pour plate technique
	Enumeration of <i>Clostridium perfringens</i>	LST EN ISO 15213-2:2024	Counting method. Pour plate technique
	Enumeration of coagulase-positive staphylococci (<i>Staphylococcus aureus</i> and other species), 37 °C	LST EN ISO 6888-1:2021, LST EN ISO 6888-1:2021/A1:2023	Counting method. Pour plate technique
	Enumeration of yeasts Enumeration of molds	LST ISO 21527-2:2008	Counting method. Principle of surface inoculation

Materials or products tested	Component, parameter or characteristic to be tested	Reference number of the document specifying test methods, clause (if relevant)	Techniques, methods and/or equipment used (where appropriate)
Drinking water	Enumeration of presumptive <i>Bacillus cereus</i> , 30 °C	LST EN ISO 7932:2005 except LST EN ISO 7932:2005/A1:2020	Counting method. Principle of surface inoculation
	Enumeration of <i>Listeria monocytogenes</i>	LST EN ISO 11290-2:2017	Counting method. Principle of surface inoculation
	Sulfite-reducing <i>Clostridium spp.</i> count	LST EN ISO 15213-1:2023	Counting method. Pour plate technique
	Enumeration of culturable microorganisms	LST EN ISO 6222:2001	Counting method. Pour plate technique
	Enumeration of <i>Pseudomonas aeruginosa</i>	LST EN ISO 16266:2008	Counting method. Principle of membrane filtration
	Enumeration of intestinal <i>Enterococci</i>	LST EN ISO 7899-2:2001	Counting method. Principle of membrane filtration
	Enumeration of cells and spores of <i>Clostridium perfringens</i>	LST EN ISO 14189:2016	Counting method. Principle of membrane filtration
	Spores of sulfite-reducing anaerobes (clostridia) count	LST EN 26461-2:2001	Counting method. Principle of membrane filtration
	Enumeration of <i>Escherichia coli</i> Enumeration of Coliforms	LST EN ISO 9308-1:2014 LST EN ISO 9308-1:2014/A1:2017	Counting method. Principle of membrane filtration
Dialysis water	Enumeration of microorganisms, 22 °C, 7 days	LST EN ISO 23500-1:2024, TGEA, PP	Counting method. Pour plate technique, TGEA

Note. In case of any discrepancies, ambiguities or disputes regarding the subject matter content between the English and Lithuanian versions of the document, the Lithuanian version shall prevail.

The accreditation certificate is signed with a qualified electronic signature as an attachment to the decision of the National accreditation body, by which it was approved