

Auto-immune disorders affect over 7% of the population. They are the 3rd largest cause of morbidity in industrialised countries (after cardiovascular diseases and cancers) and can affect people of any age.

Auto-immune disorders are due to malfunctions of the immune system which result in the manufacture of antibodies, and in some cases, lymphocytes, that act against the body's own cells. As a result, each and every organ can be affected. There are multiple factors accounting for the origin of auto-immune disorders with the involvement of genetic, endogenetic and environmental factors.



Diagnosis of auto-immune disorders

Diagnosis of an auto-immune disorder is based mainly on: clinical symptoms, immunological diagnosis and in some cases, genetic diagnosis. These debilitating disorders are often severe and at times life-threatening. Often affecting young subjects, hence the importance of correctly screening for auto-antibodies that are vital for diagnosis in the presence of relevant clinical signs.

The auto-immunity serology tests performed by Eurofins Biomnis can assist you in making a quick and reliable diagnosis to improve the treatment of your patients.

Key

Auto-antibodies with a high diagnostic value are shown in blue.

HLA genotyping can be a diagnostic support for several auto-immune diseases (for any additional information, please consult our website **www.eurofins-biomnis.com**).

Pathology ▼ Antibody ▼

nic disorders
Anti-nuclear antibodies Native anti-DNA Ab Anti-nucleosome Ab Anti-ENA Ab (Sm, RNP, SSA) Anti-Ibosomes Ab Anti-PCNA Ab Anti-C1q Ab Circulating immune complexes
Anti-nuclear antibodies Anti-histones Ab
Anti-nuclear antibodies Anti-SSA/Ro 60, SSB Ab Anti-SSA/Ro 52 Ab = anti-TRIM 21
Anti-nuclear antibodies Anti-RNP Ab
Anti-Jo1 Ab Anti-PL7 Ab, PL12, Ku, Mi2 Anti-SRP Ab Anti-Pm-Scl Ab Anti-OJ Ab, EJ, TIF/p155/140, MDA-5/CADM 140, NXP2, SAE, HMG-coA reductase
DFS 70
Anti-nucleolar Ab Anti-ScI70 Ab = topo isomerase Anti-fibrillarin Ab Anti-RNA polymerase III Ab
Anti-Pm-Scl Ab Anti-Ku Ab
Anti-centromere Ab (A/B) Anti-Th/To Ab Anti-NOR Ab
Anti-β2 glycoprotein 1 lgG/lgM Ab Anti-cardiolipin lgG/lgM Lupus anticoagulant Anti-phosphatidylethanolamine lgG/lgM Ab Anti-phosphatidylserine lgG Ab Anti-prothrombin lgG Ab Anti-annexin V lgG Ab
Anti-cyclic citrullinated peptide Ab (anti-CCP) Rheumatoid factors IgM Rheumatoid factors IgA

	Pathology ▼	Antibody ▼
		Liver
F	Primary biliary cirrhosis	Anti-mitochondrial type 2 Ab, pyruvate deshydrogenase complex = PDH Anti-gp210 Ab (nuclear pores) Anti-SP100 Ab (nuclear dots)
	Type 1 auto-immune hepatitis	Anti-smooth muscle Ab, type anti-actin Anti-soluble liver antigen Ab = SLA/LP
	Type 2 auto-immune hepatitis	Endoplasmic reticulum Ab = Anti-liver kidney microsomal Ab = anti-LKM 1) Anti-cytochrome P450 II D6 Ab Anti-cytosol Ab = anti-LC1
	Primary sclerosing cholangitis	Anti-neutrophil cytoplasmic antibodies, type X or p-ANCA
	6	Stomach
7	Type A gastritic anaemia (Biermer's anaemia)	Intrinsic factor Ab Anti-parietal cell antibodies Anti H+K+ ATPase Ab
		Pancreas
E .	Insulin-dependent diabetes	Pancreas Anti-islets of Langerhans Ab Anti-glutamate-decarboxylase Ab = GAD Anti-tyrosine-phosphatase Ab = IA2 Anti-ZnT8 Ab Anti-Insulin Ab
6	Insulin-dependent diabetes Auto-immune pancreatitis	Anti-islets of Langerhans Ab Anti-glutamate-decarboxylase Ab = GAD Anti-tyrosine-phosphatase Ab = IA2 Anti-ZnT8 Ab
		Anti-islets of Langerhans Ab Anti-glutamate-decarboxylase Ab = GAD Anti-tyrosine-phosphatase Ab = IA2 Anti-ZnT8 Ab Anti-Insulin Ab Anti-lactoferrin and carbonic anhydrase Ab
	Auto-immune pancreatitis	Anti-islets of Langerhans Ab Anti-glutamate-decarboxylase Ab = GAD Anti-tyrosine-phosphatase Ab = IA2 Anti-ZnT8 Ab Anti-Insulin Ab Anti-lactoferrin and carbonic anhydrase Ab Intestines
		Anti-islets of Langerhans Ab Anti-glutamate-decarboxylase Ab = GAD Anti-tyrosine-phosphatase Ab = IA2 Anti-ZnT8 Ab Anti-Insulin Ab Anti-lactoferrin and carbonic anhydrase Ab
	Auto-immune pancreatitis Coeliac disease = gluten	Anti-islets of Langerhans Ab Anti-glutamate-decarboxylase Ab = GAD Anti-tyrosine-phosphatase Ab = IA2 Anti-ZnT8 Ab Anti-Insulin Ab Anti-lactoferrin and carbonic anhydrase Ab Intestines Anti-transglutaminase IgA/IgG antibodies Anti-endomysium IgA/IgG Ab
	Auto-immune pancreatitis Coeliac disease = gluten intolerance	Anti-islets of Langerhans Ab Anti-glutamate-decarboxylase Ab = GAD Anti-tyrosine-phosphatase Ab = IA2 Anti-ZnT8 Ab Anti-Insulin Ab Anti-lactoferrin and carbonic anhydrase Ab Intestines Anti-transglutaminase IgA/IgG antibodies Anti-endomysium IgA/IgG Ab Anti-gliadin IgA/IgG Ab Anti-saccharomyces cerevisiae IgA/IgG Ab = ASCA
	Auto-immune pancreatitis Coeliac disease = gluten intolerance Crohn's disease	Anti-islets of Langerhans Ab Anti-glutamate-decarboxylase Ab = GAD Anti-tyrosine-phosphatase Ab = IA2 Anti-ZnT8 Ab Anti-Insulin Ab Anti-lactoferrin and carbonic anhydrase Ab Intestines Anti-transglutaminase IgA/IgG antibodies Anti-endomysium IgA/IgG Ab Anti-gliadin IgA/IgG Ab Anti-saccharomyces cerevisiae IgA/IgG Ab = ASCA Exocrine anti-pancreas Ab Anti-neutrophil cytoplasmic antibodies, type X or p-ANCA
	Auto-immune pancreatitis Coeliac disease = gluten intolerance Crohn's disease	Anti-islets of Langerhans Ab Anti-glutamate-decarboxylase Ab = GAD Anti-tyrosine-phosphatase Ab = IA2 Anti-ZnT8 Ab Anti-Insulin Ab Anti-lactoferrin and carbonic anhydrase Ab Intestines Anti-transglutaminase IgA/IgG antibodies Anti-endomysium IgA/IgG Ab Anti-gliadin IgA/IgG Ab Anti-saccharomyces cerevisiae IgA/IgG Ab = ASCA Exocrine anti-pancreas Ab Anti-neutrophil cytoplasmic antibodies, type X or

Anti-thyroglobulin Ab = Tg

myxoedema

Pathology ▼	Antibody ▼
	Kidney/Lung
Goodpasture syndrome	Anti-glomerular basement membrane Ab (not to be confused with anti-epidermal basement membrane and intracellular substance antibodies)
Granulomatosis with polyangitis (also known as Wegener's granulomatosis)	Anti-neutrophil cytoplasmic Ab (or ANCA), type anti-proteinase 3 = PR3
Eosinophilic granulomatous, microscopic polyangiitis with polyangiitis (formerly Churg and Strauss diseases), extra-capillary glomerulonephritis, etc.	Anti-neutrophil cytoplasmic antibody (or ANCA), type anti-myeloperoxydase = MPO
Idiopathic extra-membranous	Anti-PLA2R Ab

Adrena	nl/Ovaries/Testicles
Addison's disease	Anti-adrenal gland antibodies Anti-21 hydroxylase Ab
Early menopause Hypofertility	Anti-ovarian Ab Anti-testicular Ab

A CONTRACTOR OF THE PARTY OF TH	Nervous system
Myasthenia	Acetylcholine anti-receptor antibodies Anti-MuSK Ab Anti-striated muscles Ab Anti-titin Ab
Lambert-Eaton myasthenic syndrome	Anti-Voltage Gated Calcium Channel (VGCC) antibodies
Isaacs' syndrome	Anti-potassium channel Ab (VGKC), Lgi1 and Caspr2
Stiff-man syndrome	Anti-glutamate-decarboxylase Ab = GAD Anti-amphiphysin Ab
Peripheral neuropathies	Anti-myelin Ab Anti-MAG Ab Anti-gangliosides IgG/IgM Ab (GM1, GM2, GD1a, GD1b, GD3, GT1b, GQ1b, etc.)
Devic's disease	Antioptic neuromyelitis/NMO Ab = aquaporin 4
Spectrum of MOG antibody- associated disorders	Anti-MOG Ab
Paraneoplastic syndromes	Anti-neurones Ab: Hu, Ri, Yo, amphiphysin, CV2, Ma2, SOX1, zic4, Tr
Limbic encephalitis	anti-NMDAr Ab, AMPAr, VGKC, GAD, neurones

Pathology ▼	Antibody ▼

(
1	Ô	Skin
	Pemphigus	Anti-intracellular substance antibodies = desmosomes Anti-desmoglein Ab 1 and 3
	Bullous pemphigoid Gestational pemphigoid	Anti-epidermal basement membrane and intracellular substance antibodies (not to be confused with anti-glomerular basement membrane Ab) Anti-BP 180 Ab Anti-BP 230 Ab
	Epidermolysis bullosa	Anti-dermal basement membrane antibodies Anti-collagen VII Ab
	Paraneoplastic pemphigus	Anti-envoplakin antibodies
	Herpetiform dermatitis	Anti-transglutaminase IgA/IgG antibodies Anti-endomysium IgA/IgG Ab
	Subacute cutaneous lunus	Anti-SSA/Ro 60 Ah

Parathyroidism Parath	
Hypoparathyroidism	Anti-calcium sensing receptor (CaSR) Ab
	-

	Eye	
Retinopathies	Anti-recoverin Ab	

	ear/cartilage
Meniere's disease	Type 2 anti-collagen Ab
Atrophic polychondritis	Type 2 anti-collagen Ab

To find out more please contact...


Eurofins Biomnis
International Division
E-mail: international@biomnis.eurofinseu.com
www.eurofins-biomnis.com, specialities > Immunology



Biomnis

Eurofins Biomnis 17/19 avenue Tony Garnier BP 7322 – 69357 LYON Cedex 07 – FRANCE www.eurofins-biomnis.com





Autoantibodies in the diagnosis and follow-up of auto-immune diseases

